

Travel-Lite Premium Transporter 727



Feather release swingback arms.



No tools required for adjustable height back.

- High strength lightweight frame with wrap around, swing back arms.
- Multi-angle caster stem and height adjustable front caster forks.
- Padded breathable lightweight nylon upholstery, padded arms and calf pads.
- 8" X 1" caster wheels with lightweight polyurethane tires all around.
- Attendant-friendly rear locking loop brakes.
- Cam action swing away detachable footrests or legrests.
- Adjustable length foot plates with heel loops.
- High impact composite side panels, and foot plates.
- Weight limit: 16" & 18" - 260 lbs., 20" - 310 lbs.
- HCPCS Code: K1038

MODEL	SEAT			BACK	FRONT RIGGINGS	OVERALL					WEIGHT
	16"	18"	16"	15"-17"	17"- 20"	21"	34"-36"	41"	25"	10"	(w/o front riggings)
727N	16"	18"	16"	15"-17"	17"- 20"	21"	34"-36"	41"	25"	10"	19 lbs.
727NE	16"	18"	16"	15"-17"	17"- 24"	21"	34"-36"	44 1/2"	25"	10"	19 lbs.
727	18"	18"	16"	15"-17"	17"- 20"	23"	34"-36"	41"	25"	10"	20 lbs.
727E	18"	18"	16"	15"-17"	17"- 24"	23"	34"-36"	44 1/2"	25"	10"	20 lbs.
727W	20"	18"	16"	15"-17"	17"- 20"	25"	34"-36"	41"	25"	10"	21 lbs.
727WE	20"	18"	16"	15"-17"	17"- 24"	25"	34"-36"	44 1/2"	25"	10"	21 lbs.

E ELEVATING LEGRESTS N NARROW SEAT WIDTH (16") W WIDE SEAT WIDTH (20")