

Century™ LTC T8000

Long Term Care Bed

COLOR: ANTIQUE WHITE SEMI GLOSS FRAME WITH
WALNUT GRAIN HEAD & FOOT BOARDS



The Century™ Deluxe LTC Bed is specially engineered and designed for Alzheimer residents or others with chronic impairments which will reduce there degree of independence in activities of daily living. Top manufacturing technology with the highest quality components make these beds the best choice for long term care. The T8000 is a fully electric adjustable bed. The bed height, head and foot angles functions are standard along with trendelenburg and reverse trendelenburg functions that are all motor controlled with a stylish hand pendent control. Also available in Wide to meet the needs of larger patients.

Tuffcare®
HIGHER STANDARD FOR HEALTH CARE™

Century™ LTC T8000

Long Term Care Bed

COLOR: ANTIQUE WHITE SEMI GLOSS FRAME WITH WALNUT GRAIN HEAD & FOOT BOARDS

STANDARD FEATURES

- * Welded rectangular steel tubing construction with reinforced corner plates provide maximum stability.
- * Rigid high gauge steel mesh bed frame surface that provides maximum support with no mattress sagging.
- * Pre-assembled one piece grid design requires minimum setup.
- * Premier DC motor provides quiet and smooth operation, synchronized control system keeps the bed level for resident comfort and positioning.
- * Standard five-function hand pendant allows the bed to be positioned in Trendelenburg, reverse Trendelenburg and chair positions with a simple push of the button.
- * 9" clearance under recessed lower frame for easy housekeeping and accessibility for all patient handling.
- * Four 4" swivel, central locking casters conveniently located under the footboard, engage and disengage all four wheel locks with one simple step that provides increase security.
- * Standard under bed night LED light provides an added measure of security at night
- * Accommodates weights up to 600lbs.

OPTIONAL ACCESSORIES

- * Full or half rails
- * High density foam mattress
- * Alternating air pad system
- * Deluxe overbed table
- * Deluxe bedside cabinet
- * Gel overlay mattress
- * Deluxe convoluted top foam mattress
- * Alternating & low air loss mattress system
- * Patient Lift
- * Trapeze bar

SPECIFICATIONS

| MODEL | HI/LO OPERATION | HEAD/FOOT OPERATION | FABRIC HEIGHT | OVERALL WIDTH | OVERALL LENGTH | OVERALL HEIGHT | MATTRESS SIZE | BED WEIGHT | RAIL WEIGHT |
|-------|-----------------|---------------------|---------------|---------------|----------------|----------------|---------------|------------|-------------|
| T8036 | ELECTRIC | ELECTRIC | 9.75" to 31" | 36" | 87" | 19" to 46" | 36" X 80" | 205lbs | 19lbs |
| T8042 | ELECTRIC | ELECTRIC | 9.75" to 31" | 42" | 87" | 19" to 46" | 42" X 80" | 237lbs | 19lbs |



5 function control hand pendant. Also available with optional under bed LED light switch.



Optional secondary control available on side panels for patient control or on the footboard for staff control.



Color coded wheel locking pedals engage and disengage all 4 wheels with one simple step.