

Compact Featherlight™ 757

High Strength Lightweight

SEAT WIDTH AVAILABLE: 16", 18", 20"

UPHOLSTERY AVAILABLE: BLACK



Tuffcare's feather light 757, only 27 lbs. It is the best choice for person with a active life style. It is engineered to provide maximum durability and comfort for individuals who seek a modern look and great performance in a up scale wheelchair. With its adjustable front casters and rear wheel axle positions; this wheelchair offers wide selections of seat heights and angles. Come standard with quick release axle on rear wheels for transporting and storage.

Tuffcare®
HIGHER STANDARD FOR HEALTH CARE™

Compact Featherlight™ 757

High Strength Lightweight

SEAT WIDTH AVAILABLE: 16", 18", 20"

UPHOLSTERY AVAILABLE: BLACK

STANDARD FEATURES

- * High strength lightweight chrome plated frame with wrap around, swing back arms
- * Multi rear axle positions, multi-angle caster stem, & height adjustable front caster fork
- * Padded breathable nylon upholstery, padded arm pads & calf pads
- * 6" x 1" caster wheels with lightweight polyurethane tires
- * 24" x 1" rear light aluminum spoke wheels with quick release axles & polyurethane tires
- * Adjustable toggle wheel locks & rear antitippers
- * Cam action swing away detachable footrests or legrests
- * High impact composite foot plates with heel loops
- * Weight limit: 18"- 250 lbs, 20"- 280 lbs

OPTIONAL ACCESSORIES

- * Anti-theft bar
- * Wheelchair tray
- * Seat belt
- * Projection handrim
- * Anti-theft lock for front riggings
- * Pneumatic tires
- * I.V. rod attachment
- * Anti-folding device
- * Wheel lock extension tubes
- * Leg straps
- * Oxygen tank holder
- * Pneumatic tires with airless inserts

SPECIFICATIONS

MODEL	SEAT			BACK	FRONT RIGGINGS	OVERALL					WEIGHT (w/o Front Riggings)
	16"	18"-19"	16"			16"-18"	16"-21"	23"	34"-37"	43"	
757N	16"	18"-19"	16"	16"-18"	16"-21"	23"	34"-37"	43"	31"	11"	26 lbs.
757NE	16"	18"-19"	16"	16"-18"	16"-21"	23"	34"-37"	45"	31"	11"	26 lbs.
757	18"	18"-19"	16"	16"-18"	16"-21"	25"	34"-37"	43"	31"	11"	27 lbs.
757E	18"	18"-19"	16"	16"-18"	16"-21"	25"	34"-37"	45"	31"	11"	27 lbs.
757W	20"	18"-19"	16"	16"-18"	16"-21"	27"	34"-37"	43"	31"	11"	28 lbs.
757WE	20"	18"-19"	16"	16"-18"	16"-21"	27"	34"-37"	45"	31"	11"	28 lbs.

E ELEVATING LEGRESTS N NARROW 16" WIDTH



Multi position rear axle easily lowers the seat height or moves rear wheel forward or backward for better control.



High performance caster housing with stylish caster assembly and high race bearings.



Back height is adjustable between 16", 17" and 18" to best suit the patient's preference.