

# Rhino™ Lifter P420/425

Electric Lifter

FRAME COLOR IN: ANTIQUE WHITE



P425 With Full electric Adjust Base

*RHINO™ Patient Lifter is designed to help reduce the possibility of a caregiver injury and to ensure dignity in resident handling. Heavy gauge carbon steel construction with powder coated frame. State of the art electronic system, includes emergency lowering button. Accommodates all types of slings.*

**Tuffcare®**  
HIGHER STANDARD FOR HEALTH CARE™

# Rhino™ Lifter

## P420/425

Electric Lifter

FRAME COLOR IN: ANTIQUE WHITE

### STANDARD FEATURES

- \* Deluxe patient lift with powder coated carbon steel constructed frame.
- \* Low-profile base with 4 1/2" clearance allows the lift to fit under most long-term care beds.
- \* Equipped with twin 3" casters in the front and 4" locking casters in the rear.
- \* State of the art electronic system, which includes an emergency lowering button.
- \* Quick remove battery pack easily charges overnight on wall-mounted charging station.
- \* Full travel height for sling accommodates transfer to higher surfaces.
- \* New foot pedal design opens and closes the base legs with ease and eliminates annoying adjustments with handle bar.
- \* Available with full electric power base or a manual foot pedal base.
- \* Accommodates all type of slings.
- \* Accommodates weights up to 450lbs.
- \* HCPCS Code: E0635



MODEL	SLING HEIGHT	BASE CLEARANCE	BASE WIDTH	BASE LENGTH	FORK HEIGHT	WEIGHT W/O SLING	PACKAGING
P420	0" - 35"	4 1/2"	26" - 38"	47"	26" - 72"	105 lbs.	1 / CASE
P425	0" - 35"	4 1/2"	26" - 38"	47"	26" - 72"	117 lbs.	1 / CASE

### External Charging Pack PCHW1

- Designed for use as a wall charger.
- Fully charges battery pack in 4 hours.
- Automatic shut off when fully charged.
- Great for use as a spare unit or as an immediate replacement.



Emergency Operation



Removeable Battery Pack



Emergency Lowering Knob