

# Rhino™ Patient Lift P350

Hydraulic Floor Lift

FRAME COLOR IN: ANTIQUE WHITE



*RHINO™ patient Lift is designed to help reduce the possibility of a caregiver injury and to ensure dignity in resident handling. Heavy gauge carbon steel construction with powder coated frame. High performance hydraulics raise or lower individuals gradually and safely. Six point swivel cradle bar that can easily accommodate all type of sling styles and positions.*

**Tuffcare®**  
HIGHER STANDARD FOR HEALTH CARE™

# Rhino™ Patient Lift P350

Hydraulic Floor Lift


FRAME COLOR IN: ANTIQUE WHITE

## STANDARD FEATURES

- \* Deluxe manual hydraulic floor style patient lifter with powder coated carbon steel constructed frame.
- \* Low-profile base with 4 1/2" clearance allows the lift access to fit under most long-term care beds.
- \* Equipped with dual 3" front caster wheels and 4" locking rear caster wheels.
- \* High performance hydraulic cylinder with double o-ring seals provides long lasting quality and durability.
- \* Full travel height positioning for sling will accommodate patient transfer to higher areas.
- \* Adjustable base with foot pedal design allows to open and close the base legs easily with minimum effort.
- \* Six point hook cradle bar accommodates all types of slings.
- \* Accommodates weights up to 450lbs.
- \* HCPCS Code: E0630

## SPECIFICATIONS

MODEL	SLING HEIGHT	BASE CLEARANCE	BASE WIDTH	BASE LENGTH	FORK HEIGHT	WEIGHT W/O SLING	PACKAGING
P350	0" - 35"	4 1/2"	26" - 38"	47"	26" - 72"	110 lbs.	1 / CASE



Front 3" dual casters with 4 1/2" clearance.

Rear 4" casters with easy step on/off wheel lock plates.



Adjustable base with foot pedal design allows the caregiver to open and close the base legs with minimum effort for transferring and positioning safely.



Quality hydraulic manual cylinder pump with a double o-ring seal for long lasting durability which allows to raise and lower individuals safely.