



Front dual 5" casters with sealed bearing insure long ability.



Rear wheeled roller glide provides effortless sliding

- •Anodized 1" aluminum frame with dual reinforced steel cross brace.
- •Two button release mechanism.
- •Foldable for easy transporting or storage.
- •Soft, comfortable hypalon handgrips.
- Equipped with anti-rattle silencers for quiet operation.
- •Front caster wheels and rear wheeled glide for easy walking.
- •Height adjusts in 1" increments.
- Weight limit: 700lbs.HCPCS Code: E0149

MODEL	OVERALL WIDTH	HANDGRIP WIDTH	HEIGHT	OPEN DEPTH	FOLDED DEPTH	FRONT WHEEL	WEIGHT	PACKAGING PER / CASE
W500	35"	27"	34" - 38"	15"	4 1/2"	DOUBLE	16.2 lbs.	4 PER CASE