

Bariatric Commode M500



Bariatric



Drop arm mechanism is easy to operate.

- Welded polyester coated steel frame construction.
- Extra wide snap-on seat with cover.
- Hygienic seat opening.
- Plastic bucket with cover and splash guard
- Height adjusts in 1" increments.
- Blow molded armrests for extra comfort.
- Use at bedside, act as toilet safety frame or raised toilet seat.
- Weight limit: 700 lbs.
- HCPCS Code: E0168

MODEL	SEAT HEIGHT	SEAT DEPTH	ARM HEIGHT	OVERALL WIDTH	OVERALL DEPTH	WIDTH BETWEEN ARMS	WEIGHT	PACKAGING PER / CASE
M500	20" - 24"	18"	26"-31"	31"	21"	29"	22 lbs.	1 PER CASE